



Respiratory Therapy

Mini-Plus

Comfortable Compressor Nebulizer

The Mini-Plus is a compressor nebulizer designed for easy portability, especially for people on the move. Chronic respiratory problems are unpredictable and it is important for people to have their treatment devices close at hand. The compact size, light weight and convenient carrying handle of the Mini-Plus facilitate portability and allow each person to keep their Mini-Plus accessible at all times.





Features:

- It's easy to carry and ready to go right out of the box.
- Updated sleek design makes it easy to pack and store.
- Integrated handle makes it convenient for travel and mobility.
- Motor-piston drive for reliable, consistent performance.

Accessories Ordering Information:



Mouthpiece set : Anti-drool "T" mouthpiece with 6" reservoir tube and 7' supply tube.

Mouthpiece set and 2 masks : Anti-drool "T" mouthpiece with 6" reservoir tube, 7' supply tube and one each of adult & pediatric mask



Delivers aerosolized particles of 0.5-5 microns.
 Dispenses 3 cc's of medication in 10 minutes or less.
 Reliable output for efficient treatment.

Convenient Tube Holder design.

Flow outlet fitting, air filter and on/off switch are positioned around the periphery on the front panel and recessed for extra protection.

One button operation for easy set up within seconds.

Specifications:

Electrical Rating	120VAC, 60Hz, 1.2A	230VAC, 50Hz, 0.6A
(Note: Refer to the rating label on the product)	230VAC, 60Hz, 0.6A	
Maximum Compressor Pressure	≥ 35 psi	≥ 35 psi
Operation Flow rate	≥ 6.0 lpm	≥ 5.5 lpm
Nebulization Rate	0.2 ~ 0.3 ml / min	
Particle Size	0.5 ~ 5 microns	
Dimensions (W×D×H)	14.0 × 18.8 × 10.0 cm / 5.5"× 7.4"× 3.9"	
Weight	1.6 kg / 3.5 lb	



APEX MEDICAL CORP.

No. 9, Min Sheng St., Tu-Cheng,
 New Taipei City, 23679, Taiwan
 TEL:+886 2 2268 3100
 FAX:+886 2 2268 6525
 marketing@apexmedical.com.tw

APEX MEDICAL S.L.

Elcano 9, 6ª planta
 48008 Bilbao. Vizcaya. Spain
 TEL:+34 94 4706408
 FAX:+34 94 4706409
 info@apexmedical.es

www.apexbrand.com

Print-2012/08-500/All rights reserved
 Specifications may change without notice.